


Colour-marks



BioQuick® Brackets – FACE Evolution System

Maxillary					Slot .018"		Slot .022"	
Tooth	Torque	Angulation	In/Out	Rotation	Order-no.		Order-no.	
					Right	Left	Right	Left
<b>1</b> Centrals	+12°	+5°	0,8	-	739-0103	738-0103	739-0101	738-0101
<b>2</b> Laterals	+8°	+9°	1,3	-	739-0203	738-0203	739-0201	738-0201
<b>3</b> Cuspids Cuspids + hook	+3°	+8°	0,8	-	739-0303	738-0303	739-0301	738-0301
	+3°	+8°	0,8	-	739-0313	738-0313	739-0311	738-0311
<b>4</b> Bicuspids Bicuspids + hook	-7°	0°	0,9	-	739-0503	738-0503	739-0501	738-0501
	-7°	0°	0,9	-	739-0513	738-0513	739-0511	738-0511
<b>5</b> Bicuspids Bicuspids + hook	-7°	0°	0,9	-	739-0503	738-0503	739-0501	738-0501
	-7°	0°	0,9	-	739-0513	738-0513	739-0511	738-0511
<b>6</b> 1. Molar	-30°	0°	-	10	739-0704	738-0704	739-0702	738-0702
<b>7</b> 2. Molar	-30°	0°	-	6	739-0804	738-0804	739-0802	738-0802

Mandibular					Slot .018"		Slot .022"	
Tooth	Torque	Angulation	In/Out	Rotation	Order-no.		Order-no.	
					Right	Left	Right	Left
<b>1</b> Centrals	-1°	0°	1,4	-	738-1303	738-1303	738-1301	738-1301
<b>2</b> Laterals	-1°	0°	1,4	-	738-1303	738-1303	738-1301	738-1301
<b>3</b> Cuspids Cuspids + hook	-6°	+2°	1,2	-	739-1403	738-1403	739-1401	738-1401
	-6°	+2°	1,2	-	739-1413	738-1413	739-1411	738-1411
<b>4</b> Bicuspids Bicuspids + hook	-17°	0°	1,2	-	739-1503	738-1503	739-1501	738-1501
	-17°	0°	1,2	-	739-1513	738-1513	739-1511	738-1511
<b>5</b> Bicuspids Bicuspids + hook	-22°	0°	1,2	-	739-1603	738-1603	739-1601	738-1601
	-22°	0°	1,2	-	739-1613	738-1613	739-1611	738-1611
<b>6</b> 1. Molar	-30°	0°	-	-	739B1704	738B1704	739B1702	738B1702
<b>7</b> 2. Molar	-30°	0°	-	-	739-1804	738-1804	739-1802	738-1802

Fallsätze | Cases c € 0297

Variation	Slot .018"			Slot .022"		
	1	5	10	1	5	10
	706-1200	706-1201	706-1202	706-1209	706-1210	706-1211
3	706-1203	706-1204	706-1205	706-1212	706-1213	706-1214
3-5	706-1206	706-1207	706-1208	706-1215	706-1216	706-1217

Optionally more torque on mandibular Centrals and Laterals

Mandibular					Slot .018"		Slot .022"	
Tooth	Torque	Angulation	In/Out	Rotation	Order-no.		Order-no.	
					Right	Left	Right	Left
<b>1+2</b> Centrals / Laterals	-6° *	0°	1,4	-	738-1203	738-1203	738-1201	738-1201

\* Bonding brackets upside down will convert from -6° Torque to +6° Torque. But keep in mind: the clip has to be opened to gingival.  
Please indicate if you prefer this option.