


Colour-marks



Quicklear® Brackets – FACE Evolution System

Maxillary <span style="float:right">C€ 0297</span>					Slot .018"		Slot .022"	
Tooth	Torque	Angulation	In/Out	Rotation	Order-no.		Order-no.	
					Right	Left	Right	Left
<b>1</b> Centrals	+12°	+5°	0,8	-	C739-0103	C738-0103	C739-0101	C738-0101
<b>2</b> Laterals	+8°	+9°	1,3	-	C739-0203	C738-0203	C739-0201	C738-0201
<b>3</b> Cuspids + hook	+3°	+8°	0,8	-	C739-0313	C738-0313	C739-0311	C738-0311
<b>4</b> Bicuspids	-7°	0°	0,9	-	C739-0503	C738-0503	C739-0501	C738-0501
	Bicuspids + hook	-7°	0°	0,9	-	C739-0513	C738-0513	C739-0511
<b>5</b> Bicuspids	-7°	0°	0,9	-	C739-0503	C738-0503	C739-0501	C738-0501
	Bicuspids + hook	-7°	0°	0,9	-	C739-0513	C738-0513	C739-0511
<b>6</b> 1. Molar	-30°	0°	-	10	739-0704	738-0704	739-0702	738-0702
<b>7</b> 2. Molar	-30°	0°	-	6	739-0804	738-0804	739-0802	738-0802

Mandibular <span style="float:right">C€ 0297</span>					Slot .018"		Slot .022"	
Tooth	Torque	Angulation	In/Out	Rotation	Order-no.		Order-no.	
					Right	Left	Right	Left
<b>1</b> Centrals	-1°	0°	1,4	-	C738-1303	C738-1303	C738-1301	C738-1301
<b>2</b> Laterals	-1°	0°	1,4	-	C738-1303	C738-1303	C738-1301	C738-1301
<b>3</b> Cuspids + hook	-6°	+2°	1,2	-	C739-1413	C738-1413	C739-1411	C738-1411
<b>4</b> Bicuspids	-17°	0°	1,2	-	739-1503	738-1503	739-1501	738-1501
	Bicuspids + hook	-17°	0°	1,2	-	739-1513	738-1513	739-1511
<b>5</b> Bicuspids	-22°	0°	1,2	-	739-1603	738-1603	739-1601	738-1601
	Bicuspids + hook	-22°	0°	1,2	-	739-1613	738-1613	739-1611
<b>6</b> 1. Molar	-30°	0°	-	-	739B1704	738B1704	739B1702	738B1702
<b>7</b> 2. Molar	-30°	0°	-	-	739-1804	738-1804	739-1802	738-1802

Fallsätze | Cases C€ 0297

Variation	Slot .018"			Slot .022"		
	1	3	6	1	3	6
3	C706-1203	C706-1204	C706-1205	C706-1212	C706-1213	C706-1214
3-5	C706-1206	C706-1207	C706-1208	C706-1215	C706-1216	C706-1217

Optionally more torque on mandibular Centrals and Laterals

Mandibular <span style="float:right">C€ 0297</span>					Slot .018"		Slot .022"	
Tooth	Torque	Angulation	In/Out	Rotation	Order-no.		Order-no.	
					Right	Left	Right	Left
<b>1+2</b> Centrals / Laterals	-6° *	0°	1,4	-	C738-1203	C738-1203	C738-1201	C738-1201

\* Bonding brackets upside down will convert from -6° Torque to +6° Torque. But keep in mind: the clip has to be opened to gingival.  
Please indicate if you prefer this option.